

Glass body

HEAVY DUTY pH ELECTRODE

Model : PE-02

ISO-9001, CE, IEC1010



← Glass body



Lutron

LUTRON ELECTRONIC

The Art of Measurement

Glass body, Professional

HEAVY DUTY pH ELECTRODE

Model : PE-02

Features	Professional, heavy duty and wide Temp. range pH electrode, available for the professional, Laboratory, high quality and field pH measurement.	
Measuring Range	0 to 14 pH .	
Operation Temperature	5 to 90 °C (41 to 194 °F).	
Electrode Structure	Combination type.	
Zero Potential for pH Value	7 ± 1 pH.	
Pressure rating	6 bar max.	
Body Material	Glass	
Connector	BNC	
Weight	75 g/0.17 LB.	
Dimensions	Body length	155 mm
	Body Dimension	15 mm
	Cable length	1 meter
Applications	Water conditioning, aquariums, beverage, fish hatcheries, food processing, photography, laboratory, paper industry, plating industry, quality control, school & college, water conditioning.	